

---

# **SARS-CoV-2 Rapid Antigen Test**

*Primary ad secondary influencers on assay performance*

*11-Dec-2020*

# Objectives of this presentation

- The main objective is to summarize key publications that deal with real world performance of the SARS-CoV-2 Rapid Antigen Test
- Secondary objective to showcase factors that influence assay performance
- This presentation will be updated regularly
- Literature search criteria and outcome will be listed
- Publications with results that include comparisons with an CE /EUA approved PCR and the corresponding Ct values will be summarized in this presentation

Roche SARS-CoV-2 Rapid Antigen Test = STANDARD Q-COVID-19 Ag Test


# Search strategy; 30 – Nov - 2020

Set#	Searched for	Results
S1	(Ti,Ab(COVID-19 OR "COVID-19" OR COVID19 OR SARS-CoV-2 OR SARSCoV2 OR SARS-CoV2 OR "SARS-CoV-2")) OR (MJEMB.EXACT.EXPLODE("severe acute respiratory syndrome")) OR MJEMB.EXACT.EXPLODE("Coronaviridae") OR MJMESH.EXACT.EXPLODE("Coronaviridae") OR (MJMESH.EXACT.EXPLODE("Severe Acute Respiratory Syndrome"))	179395
S2	emb("coronavirus disease 2019 +")	66193
S3	(((novel NEAR/5 corona NEAR/5 virus) OR (2019 NEAR/2 nCoV) OR ((2019 or novel) NEAR/2 coronavirus*) or "2019-nCoV" or "COVID-19" or (COVID PRE/0 19) or (corona NEAR/5 virus NEAR/5 2019) or (SARS pre/0 CoV pre/0 2) or "SARS-CoV-2"))	170265
S4	S3 OR S2 OR S1	194253
S5	("STANDARD Q COVID-19 Ag")	4
S6	(rapid n/5 antigen* n/5 (test* or assay*))	5886
S7	((S5 or S6) and S4)	70°
S8	(EMB.EXACT.EXPLODE("point of care testing")) OR (MESH.EXACT.EXPLODE("Point-of-Care Testing")) OR (poc or point n/2 care)	90332*
S9	(s4 and s8)	888°
S14	(ti,ab,su,emb,mesh(clinical n/2 perform*)) OR (ti,ab,su,emb,mesh(accuracy* OR sensitiv* OR specific* OR validation* OR concordance* OR "positive agreement" OR "positive percent agreement" OR "negative agreement" OR "negative percent agreement" OR evaluat* OR performance* OR "clinical performances"))	27646286*
S15	(s7 and s14) (ausgeliefert)	48°
S16	(s9 and s14)	471°
S17	((s9 and s14)) and (pd(20190101-20211231))	460°
S18	(s17 not s15) => zusätzliche Publikationen, gefunden mit PoC (Point of Care)	444°

\* Duplicates are removed from the search, but included in the result count.  
 ° Duplicates are removed from the search and from the result count.

## Databases:

- BIOSIS Previews®
- Derwent Drug File
- Embase®
- MEDLINE®



# Factors impacting on Performance and Test Results of Rapid Antigen Tests

# Coronaviruses

## *Virion morphology and structural proteins*

### Large enveloped RNA viruses (80–120 nm)<sup>1-3</sup>

#### Lipid bilayer



Membrane protein (M)

Envelope protein (E)

Spike protein (S)

#### Nucleocapsid


Multiple copies of the **nucleocapsid protein (N)** bound to the RNA genome



1. Masters PS (2006). Advances in Virus Research. Academic Press. 66: 193–292; 2. Su, S et al. (2016). Trends in Microbiology. 24 (6): 490–502; 3. Paules CI et al. (2020). JAMA. 2020;323(8):707–708

# Summary: Factors Impacting on Performance and Test Results of Rapid Antigen Tests

## Primary influencer:



**Viral load** of the sample, and the **viral load distribution** in the investigated cohort represented by Cycle threshold (Ct) of the PCR




**Analytical test performance of the assay:** sensitivity & specificity

## Secondary influencer:



**Days post symptom onset (DPSO)** of sampling



**Pretest probability or prevalence setting of test**




**Sampling method, e.g.**

- Swabs
- Tubes
- Buffer, Viral Transport Media



**Sample Type**

- Naso-/Oropharyngeal
- Nasal
- Saliva



**Workflow**

- Point of Care setting
- Laboratory
- Storage

1. Wölfel R, Corman VM, Guggemos W, et al. Virological assessment of hospitalized patients with COVID-2019. Nature. 2020;581(7809):465-469. doi:10.1038/s41586-020-2196-x

2. Krueger et al, <https://www.medrxiv.org/content/10.1101/2020.10.01.20203836v1>; 3. Van Beek, J et al: <https://doi.org/10.1101/2020.10.13.20211524> ; 4. Lee R. et al. Performance of Saliva, Oropharyngeal Swabs, and Nasal Swabs for SARS-CoV-2 Molecular Detection: A Systematic Review and Meta-analysis medRxiv 2020.11.12.20230748; doi: <https://doi.org/10.1101/2020.11.12.20230748>

# Influencers of Test Performance



**Sampling type/  
specimen source**

+



**Collection device /  
Transport media and  
volume**

+




**Time to test /  
transport /  
storage**

+



**Test type /  
target**

=



**Viral load of the sample /  
distribution in a cohort**


Days from infection to specimen collection

**Pre-analytical**

**Analytical**

1. Wölfel R, Corman VM, Guggemos W, et al. Virological assessment of hospitalized patients with COVID-2019. Nature. 2020;581(7809):465-469. doi:10.1038/s41586-020-2196-x  
2. Krueger et al, <https://www.medrxiv.org/content/10.1101/2020.10.01.20203836v1>; 3. Van Beek, J et al: <https://doi.org/10.1101/2020.10.13.20211524>; 4. Lee R. et al. Performance of Saliva, Oropharyngeal Swabs, and Nasal Swabs for SARS-CoV-2 Molecular Detection: A Systematic Review and Meta-analysis medRxiv 2020.11.12.20230748; doi: <https://doi.org/10.1101/2020.11.12.20230748>

# Clinical Sensitivity of a Rapid Test compared to PCR



\*Of note, Ct values are not directly translatable between different PCR methods; even the technical limit of detection can vary greatly among the EUA-approved PCR platforms. Thus the Ct value comparison here rather illustrates a trend and is not precise

\*\*Curve is for illustrative purposes only

WHO update webinar Sept 11, 2020

Wölfel et al 2020, <https://doi.org/10.1038/s41586-020-2196-x>


# Targets of different Rapid Ag tests




Different assays target different components of the SARS-CoV-2



Targets the **Nucleocapid**




Even with the same target, the antibodies may have **different epitopes and affinities**




# Quality of Samples for COVID-19 Testing

*Viral load differs for sample types and different disease severities*

## Time from symptom onset



## Disease severity



**Viral load on swabs decreases as symptoms resolve or disease progresses into lungs**


**Higher viral loads associated with more severe disease**

1. Wölfel R, Corman VM, Guggemos W, et al. Virological assessment of hospitalized patients with COVID-2019. Nature. 2020;581(7809):465-469. doi:10.1038/s41586-020-2196-x | 2. Magleby R, Westblade LF, Trzebucki A, et al. Impact of SARS-CoV-2 Viral Load on Risk of Intubation and Mortality Among Hospitalized Patients with Coronavirus Disease 2019 [published online ahead of print, 2020 Jun 30]. Clin Infect Dis. 2020;ciaa851. doi:10.1093/cid/ciaa851


# Quality of Samples for COVID-19 Testing

## *Viral load differs across storage conditions*

### Collection media and swab



### Storage temperature and time



**Test samples as soon as possible after collection**

**To improve detection, store samples refrigerated and/or in buffered viral transport media containing antibiotics**


## Stability of viral RNA affected by collection media and storage conditions

| 1. Kim N, Kwon A, Roh EY, et al. Effects of Storage Temperature and Media/Buffer for SARS-CoV-2 Nucleic Acid Detection [published online ahead of print, 2020 Oct 17]. Am J Clin Pathol. 2020;aqaa207. doi:10.1093/ajcp/aqaa207 | 2. Druce J, Garcia K, Tran T, Papadakis G, Birch C. Evaluation of swabs, transport media, and specimen transport conditions for optimal detection of viruses by PCR. J Clin Microbiol. 2012;50(3):1064-1065. doi:10.1128/JCM.06551-11

# Is a quantitative test (viral load) useful?



**Cycle threshold (C<sub>t</sub>):** Number of PCR cycles needed to produce a positive result



Lower C<sub>t</sub> value

=



Higher concentrations of viral RNA in the sample



No quantitative SARS-CoV-2 assays have received Emergency Use Authorization (EUA) by the Food and Drug Administration (FDA).




No international, commutable standardized reference material is currently available

# Comparing sensitivities of SARS-CoV-2 rapid antigen tests



**Sensitivities** of rapid antigen tests can only be compared:




An absolute assessment of **limits of detection** for each test, as well as a strict comparison of **relative sensitivities** is **not possible**

# Comparing Ct values

*Ct-values can vary significantly between and within methods*

## CAP survey

>700 laboratories using proficiency testing material produced from the same batch



Rhoads et al 2020 ; DOI: 10.1093/cid/ciaa1199



### Different FDA EUA methods:

Median Ct-values for varied by as much as **14 cycles**



### Different targets - one instrument:


Within a single test performed, the difference in the median Ct-values for different targets was **3.0 cycles**



### Across all labs:

Within a single gene target for a single method, up to **12.0 cycle** differences  
ORF1a detection differed by **6.0 cycles**



The background features a dark blue color scheme with various chemical structures and molecular models. A prominent feature is a large, glowing blue virus particle with a spherical head and a tail-like structure, resembling a bacteriophage or a similar virus. The overall aesthetic is scientific and technological.

# External Clinical Validation Studies of the Roche SARS-CoV-2 Rapid Antigen Test

(= *STANDARD Q-COVID-19 Ag Test*)

# FIND REPORT: Summary



## Purpose of the study

Independent evaluation of the performance of the test in different patient populations and prevalence settings, performed in three independent sites, two in Germany (Heidelberg and Berlin) and one in Brazil (Macaé, state of Rio de Janeiro). Patients included in the study were those that fulfilled the respective national suspect definition at the time of the study.

## Main results

Combined overall sensitivity was 84.97% with a specificity of 98.84%.  
The combined sensitivity for  $Ct \leq 25$  was 97.14%.

## Specifics

This study was designed according to the requirements of WHO Emergency Use Listing (EUL). The two German cohorts and the Brazilian cohort have to be viewed as one study, as neither site / country would fulfill these criteria alone. The WHO EUL of SD Biosensor is also based on the combined data (Germany & Brazil combined).

## Main Conclusions

The Roche SARS-CoV-2 Rapid Antigen Test is a reliable test providing fast answers wherever they are needed

FIND data complement the IFU data and give more information about the performance of the test in different settings.

# FIND REPORT: Patient Characteristics\*



	Germany	Brazil
<b>N, PCR + (%)</b>	1259 (3.7%)	400 (26.5%)
<b>Investigated cohort</b>	symptomatic & asymptomatic meeting national <suspect> definition	symptomatic & asymptomatic meeting national <suspect> definition
<b>Study + sample size</b>	Nasopharygeal and oropharyngeal	Nasopharyngeal
<b>Symptomatics, n (%)</b>	1039 (84.7%)	392 (98.7%)
DPSO (median (Q1-Q3))	3 (2-4)	5 (4-6)
Days < 0-3)	62.7%	21.4%
Days 4-7	30.9%	68.8%
Days 8+	6.4%	9.8%
<b>PCR Ct (median)</b>	<b>25.3</b>	<b>25.5</b>
CT > 33 (n,%)	6 (12.8%)	7 (6.6)
CT > 30 (n,%)	11 (23.4%)	19 (17.9%)
CT > 25 (n,%)	26 (55.3%)	57 (53.8%)
<b>Reference Method</b>	1. cobas 2. Abbott 3. Genesig (Primerdesign) 4. Allplex (Seegane) 5. LightMix (Tib Molbiol)	1. Lab-developed assays based on US CDC protocol, which targets 2 regions (N1+N2) of the NC gene (FDA EUA)

\*fulfilling WHO requirements on Emergency Use Listing (EUL)

[https://www.finddx.org/wp-content/uploads/2020/09/SDQ-Ag-Public-Report\\_20200918.pdf](https://www.finddx.org/wp-content/uploads/2020/09/SDQ-Ag-Public-Report_20200918.pdf)

# FIND REPORT: Assay Performance



 **Combined**

 **Germany**

 **Brazil**

<b>Sensitivity Ct ≤ 25</b>	<b>97.14%</b> (95% CI 90.1% – 99.65%)	<b>100%</b> (95% CI 84.5% – 100%)	<b>95.9%</b> (95% CI 86.3% – 95.9%)
<b>Sensitivity Ct ≤ 33</b>	<b>90.7%</b> (95% CI 84.6% – 95%)	<b>87.8%</b> (95% CI 74.5% – 94.7%)	<b>91.9%</b> (95% CI 84.9% – 95.9%)
<b>Sensitivity ≤ 7 days (85% CI)</b>	<b>87.88%</b> (95% CI 81.06% – 92.9%)	<b>80%</b> (95% CI 64.1% – 90.1%)	<b>90.7%</b> (95% CI 74.583.3 – 95.0%)
<b>Sensitivity (95% CI)</b>	<b>84.97%</b> (95% CI 78.3% – 90.23%)	<b>76.6%</b> (95% CI 62.8% – 86.4%)	<b>88.7%</b> (95% CI 81.3% – 93.4%)
<b>Specificity</b>	<b>98.94%</b> (95% CI 98.23% – 99.39%)	<b>99.3%</b> (95% CI 98.6% – 99.6%)	<b>97.6%</b> (95% CI 95.2% – 98.8%)

# FIND REPORT: Differences between the two cohorts



3,7% of the German cohorts and 26,5% of the Brazilian cohort tested positive by PCR.

84,7% of the German cohorts and 98,7% of the Brazilian cohort were symptomatic.

The median days post symptom onset (DPSO) is slightly lower in the German cohorts (3 DPSO) than in the Brazilian cohort (5 DPSO).

Different PCR reference methods were used (Ct values are not comparable as RT-PCR methods vary across sites with different genome targets, PCR instruments and reagents).

The two sites in Germany had more low viral-load samples (23,4% of Ct > 30; 12,8% Ct > 33) than the site in Brazil (17,9% Ct > 30; 6,6% Ct > 33)

For some patients in the study oropharyngeal swabs were used (not NP) which is not according the IFU.

[https://www.finddx.org/wp-content/uploads/2020/09/SDQ-Ag-Public-Report\\_20200918.pdf](https://www.finddx.org/wp-content/uploads/2020/09/SDQ-Ag-Public-Report_20200918.pdf)

# Hospital Universitaires Genève (HUG), Switzerland: Study Summary



## Purpose of the study

SARS-CoV-2 antigen rapid diagnostic test (RDT) validation for Panbio™ Covid-19 Ag Rapid Test (Abbott) and Standard Q COVID-19 Rapid Antigen Test (SD Biosensor/Roche), partly done in collaboration with the Foundation for Innovative Diagnostics (FIND), Geneva and supported by the CRIVE and The Geneva Centre for Emerging Viral Diseases

## Main results

RDT test results show highest concordance in samples with low CT values (indicating a high viral load). The overall sensitivity was 89%, for Ct values between <26 it was 90-100%. Despite more samples with lower viral load, Roche Ag Test shows better overall sensitivity and esp. for Ct values 26 – 48 (low viral load).

## Specifics

First swab was used for PCR, second for the Rapid Antigen testing. Second swabs might contain lower viral load.

This report will be completed as a full paper rapidly.

## Main Conclusions

The results show that the Standard Q (SD Biosensor/Roche), fulfil the criteria as defined by WHO with 80% sensitivity and 97% specificity , which is in line with independent validations from other studies.

# Hospital Universitaires Genève, Switzerland: Study Details



**Roche Ag Test**



**Abbott PanBio**



N, PCR + (%)	529 <b>(36%)</b>	535 <b>(23%)</b>
<b>Investigated cohort</b>	Symptoms for 0-4 days, n (%) 141, (77%) symptomatic & asymptomatic meeting national <suspect> definition	Symptoms for 0-4 days, n (%) 86, (75,4%) symptomatic & asymptomatic meeting national <suspect> definition
<b>Samples</b>	Nasopharyngeal, 1. swab for PCR, 2. swab for POC test	Nasopharyngeal, 1. swab for PCR, 2. swab for POC test
<b>Sensitivity overall</b>	<b>89.0%</b> (95% CI 83.69-93.06)	<b>85.48%</b> (95% CI 78.03-91.16%)
Symptoms for 0-4 days	90.85%	87.21%
Ct 14- 18	100%	96%
Ct 18-22	98%	100%
Ct 22-26	90%	90%
Ct 26-30	84%	64%
Ct 30-34	45%	38%
Ct 34-48	17%	0%
<b>Specificity</b>	<b>99.70%</b> (95%CI 98.36-99.99)	<b>100%</b> (95% CI 99.11-100.0)
<b>Positive Predictive Value</b>	<b>99.42%</b> (95%CI 96.00-99.92)	<b>100%</b>
<b>Negative Predictive Value</b>	<b>94.13%</b> (95%CI 91.47-96.00)	<b>95.80%</b> (93.71-97.22)
<b>Reference Method</b>	cobas, Roche	cobas, Roche

[https://www.hug.ch/sites/interhug/files/structures/laboratoire\\_de\\_virologie/documents/Centre\\_maladies\\_virales\\_infectieuses/ofsp\\_rdt\\_report\\_gcevd\\_27.10.2020.pdf](https://www.hug.ch/sites/interhug/files/structures/laboratoire_de_virologie/documents/Centre_maladies_virales_infectieuses/ofsp_rdt_report_gcevd_27.10.2020.pdf)


# Hospital Universitaires Genève: Result Details



## Rapid Test Sensitivity



Roche Ag Test with more samples with lower viral load and higher sensitivity for Ct values 26 - 48



[https://www.hug.ch/sites/interhug/files/structures/laboratoire\\_de\\_virologie/documents/Centre\\_maladies\\_virales\\_infectieuses/ofsp\\_rdt\\_report\\_gcevd\\_27.10.2020.pdf](https://www.hug.ch/sites/interhug/files/structures/laboratoire_de_virologie/documents/Centre_maladies_virales_infectieuses/ofsp_rdt_report_gcevd_27.10.2020.pdf)

# Cerutti et al., Italy: Study Summary



## Purpose of the study

This study evaluated the sensitivity, specificity, negative and positive predictive values (NPV and PPV) of the STANDARD Q COVID-19 Ag point-of-care diagnostic test (POCT) for the detection of SARS CoV-2 nucleoprotein in nasopharyngeal swab, in comparison with the gold standard RT-PCR

## Main results

The STANDARD Q COVID-19 Ag test showed an overall 70.6 % sensitivity and 100% specificity presenting with a Ct between 12.3 - 38.5. For samples with a Ct < 28 the sensitivity was 100%. Screening of asymptomatic persons without contact to a confirmed case results in lower performance.

## Specifics

A major limit of the study was that the test was assessed in suboptimal conditions using UTM samples instead of on-site NP swabs.

Ct values and categories are not comparable with other studies. 3 different PCR methods were used.

## Main Conclusions

The POC test shows good sensitivity for investigation of symptomatic patients. POCT (discrepant to PCR) negative results were found in samples with a low viral load, consistent with low viable virus and low infectiousness as confirmed by cell-culture in a subset of samples.

# Cerutti et al., Italy: Study Details



## Diagnostic Population 1

## Screening Population 2

<b>N, PCR positive (%)</b>	330 <b>(33%)</b>	
<b>N, PCR positive (%)</b>	185 <b>(56%)</b>	145 <b>(3.4%)</b>
<b>Investigated cohort</b>	185 with symptoms and signs consistent with COVID-19	145 asymptomatic travelers returning from EU high risk countries
<b>Samples</b>	Nasopharyngeal (NP), COPAN UTM; A major limit of the study was that the test was assessed in suboptimal conditions using UTM samples instead of on-site NP swabs. 13/185, 7% Ag tests were run on left-over sample stored at -20 °C.	
<b>Sensitivity</b>	<b>72.1%</b>	<b>40%</b>
<b>Sensitivity overall</b>	<b>70.6%</b>	
<ul style="list-style-type: none"> <li><b>Sensitivity at Ct &lt;28</b></li> <li>Ct 28 - 30</li> <li>Ct 30 - 35</li> <li>Ct &gt; 35</li> </ul>	<b>100%</b> 38.5% 26.7% 9.1%	
<b>Specificity, positive/total nr</b>	<b>100%</b> (81/81)	<b>100%</b> (140/140)
<b>Positive Predictive Value</b>	<b>100%</b>	<b>100%</b>
<b>Negative Predictive Value</b>	<b>73.6%</b>	<b>97.9%</b>
<b>Reference Method</b>	SeegeneAllplex (n=159), cobasRoche (n=118), DiaSorinSimplexa (n=28)	

Ct values not well comparable with other studies

UTM, viral transport media

Cerutti F, Burdino E, Milia MG, et al. Urgent need of rapid tests for SARS CoV-2 antigen detection: Evaluation of the SD-Biosensor antigen test for SARS-CoV-2 [published online ahead of print, 2020 Sep 29]. *J Clin Virol.* 2020;132:104654. doi:10.1016/j.jcv.2020.104654

# Krueger et al., Germany: Study Summary



## Purpose of the study

Evaluation of the accuracy, ease of use and limit of detection of novel, rapid, antigen-detecting point-of-care diagnostics for SARS-CoV-2.

Performance of three Ag-RDTs was compared to RT-PCR overall, according to predefined subcategories e.g. cycle threshold (CT)-value, days from symptoms onset. (Berlin, Heidelberg and Liverpool)

## Main results

There is large variability on performance of rapid antigen tests.

The Roche / SDB STANDARD Q-CoV test was the best performing, with 100% sensitivity for samples with Ct values < 25 and with 76.6% overall sensitivity.

## Specifics

For some patients in the study oropharyngeal samples swabs were used (not nasopharyngeal) which is not according the IFU.

The test was considered easy-to-use and suitable for point-of-care.

## Main Conclusions

With a sensitivity of 100% for the STANDARD Q COVID-19 Ag test in infected persons with a high viral load, it is likely to identify highly contagious individuals.

The rapid turn-around time is likely to result in more rapid isolation of cases and effective contact tracing.

# Krueger et al., Germany: Study Details



	Roche SARS-CoV-2* Rapid Ag	Bioeasy 2019-nCoV Ag	CorisRespi-Strip
<b>N, PCR positive (%)</b>	1263 <b>(3%)</b>	729 <b>(2.9%)</b> ,	425 <b>(1.9%)</b> ,
<b>Investigated cohorts</b>	84.4% symptomatics	81.2% symptomatics	68.9% symptomatics
<b>Samples</b>	Nasopharyngeal and oropharyngeal	Nasopharyngeal	Nasopharyngeal
<b>Sensitivity (95% CI)</b>	<b>76.6%</b> (62.8-86.4)	<b>66.7%</b> (41.7-84.8)	<b>50%</b> (21.5-78.5)
<ul style="list-style-type: none"> <li>• <b>Sensitivity</b></li> <li>• Ct &lt; 25, (95%CI)</li> <li>• Ct ≥ 25, (95%CI)</li> </ul>	<ul style="list-style-type: none"> <li>• <b>100%</b> (82.4-100)</li> <li>• <b>62.1%</b> (44.0-77.3)</li> </ul>	<ul style="list-style-type: none"> <li>• <b>88.9%</b> (56.5-99.4)</li> <li>• <b>33.33%</b> (9.7-70.0)</li> </ul>	<ul style="list-style-type: none"> <li>• <b>66.7%</b> (20.8-98.3)</li> <li>• <b>40%</b> (11.8-76.9)</li> </ul>
<b>Specificity (95%CI)</b>	<b>99.3%</b> (98.6-99.6)	<b>93.1</b> (91.0-94.8)	<b>95.8 (93.4-97.4)</b>
<b>Reference Method</b>	TibMolbiol, Allplex Seegene, Abbott, cobas 6800/8800, Genesig (UK)		

\*This is partially the data of the German cohort in the FIND study.

Krueger et al, <https://www.medrxiv.org/content/10.1101/2020.10.01.20203836v1>

# Van Beek et al., The Netherlands: Study Summary



## Purpose of the study

Freshly collected nasal and nasopharyngeal samples in viral transport media from people presenting to the drive through test station with a range of Ct values were tested in parallel by RT-PCR, and rapid antigen detection tests (RDT). Detection limits of 5 commercially available RDT's were determined using serial dilutions of freshly harvested SARS-CoV-2 virus stock.

## Main results

Rapid antigen tests differ greatly in their ability to detect infectious cases. The test were classified into 3 performance categories without further details  
With the most sensitive RDTs, 97.3% of potentially infectious individuals with mild symptoms would be detected, with medium quality tests 92.73% and with the low quality 75.53%.

## Specifics

Routine application of rapid antigen testing increased time-to-result at same day from 33% to 97%.  
Freshly collected nasal + nasopharyngeal samples in VTM tested by RT-PCR and RDT in parallel. In addition, some samples were also used for virus culture on Vero E6 cells.

## Main Conclusions


The use of rapid antigen tests for screening of individuals offers the potential for rapid identification of those individuals at greatest risk of spreading the infection. High quality RDTs offer hope to improve containment by more rapid isolation and contact tracing of the most infectious individuals.

# Van Beek et al., The Netherlands: Detection of culture positive (RT-PCR-confirmed) cases by rapid antigen tests depending on severity of symptoms

Rapid Antigen Assay	Mild, outpatient Median (min - max)	Hospitalised, mild Median (min - max)	Hospitalised, severe Median (min - max)	Roche & Abbott assays
<b>A</b> - Panbio™ COVID-19 Ag rapid test (Abbott), and <i>Standard Q COVID-19 Ag (SD Biosensor)</i>	<b>94.30%</b> (88.65% - 99.77%)	<b>98.68%</b> (95.79% - 99.81%)	<b>99.80%</b> (99.32% - 99.97%)	
<b>B</b> - COVID-19 Ag Respi-Strip (Coris BioConcept, and <i>GenBody COVID-19 Ag (GenBody Inc)</i>	<b>92.73%</b> (60.30% - 99.77%)	<b>97.43%</b> (86.40% - 99.81%)	<b>99.54%</b> (97.45% - 99.97%)	
<b>C</b> - Biocredit COVID-19 Ag (RapiGEN)	<b>75.53%</b> (17.55% - 99.75%)	<b>91.70%</b> (57.90% - 99.81%)	<b>98.55%</b> (88.53% - 99.97%)	

Rapid Antigen Tests Performance Comparison including virus culture testing of infectiousness

# Van Beek et al., The Netherlands : Correlation of PCR-/AG-test positive and cell-culture positive result for different rapid AG test performance assays



Distribution of viral RNA loads at time of diagnosis with RT-PCR confirmed SARS-CoV-2 infection N=1754 (of which 78 were tested by virus culture).

Van Beek, J et al: <https://doi.org/10.1101/2020.10.13.20211524>

# Corman et al., Germany: Study Summary



## Purpose of the study

7 different Ag POC tests were evaluated on recombinant nucleoprotein, cultured endemic and emerging coronaviruses, stored clinical samples with known SARS-CoV-2 viral loads (n=138), stored samples from patients with respiratory agents other than SARS-CoV-2 (n=100), as well as self-sampled swabs from healthy volunteers (n=35).

## Main results

The sensitivity range of most AgPOCT overlaps with viral load figures typically observed during the first week of symptoms, which marks the infectious period in the majority of patients.

All tests x-react with SARS-CoV

## Specifics

Specimens were stored in universal transport medium (Copan UTM™) at -20°C. They used stored swabs obtained in universal transport medium (Copan UTM™) or without any medium (dry swabs).

Healthy volunteers (for specificity testing) conducted self-testing. They refer to Krueger that show equivalence of specimen material.

## Main Conclusions

In hospitalized patients at the end of their clinical course, negative AgPOCT results may provide an additional criterion to safely discharge patients. Novel public health concepts suggest decisions to isolate or maintain isolation that are based on infectivity testing rather than infection screening.

# Corman et al., Germany: Study Details



## Roche Rapid Ag Test



## Abbott PanBio

N, PCR + (%)	N=529 (archive specimen)	N=535 (archive specimen)
<b>Investigated cohort</b>	symptomatic & asymptomatic meeting national <suspect> definition	symptomatic & asymptomatic meeting national <suspect> definition
<b>Samples</b>	Nasopharyngeal, swabs, dry swabs Specimens were stored at -20°C in phosphate-buffered saline (PBS) or universal transport medium (Copan UTM™) at -20°C. For specificity: self-testing	
<b>Sensitivity overall</b>	6.78 x10 <sup>6</sup> copies/swab LoD, 95% mean hit rate 4.4 PFU of virus per test	6.55 x10 <sup>6</sup> copies/swab 4.4 PFU of virus per test
<b>Specificity Cumulative Specificity</b>	<b>97.12% n= 35 98.53%</b>	<b>100% n=35 99.26%</b>
<b>Positive Predictive Value</b>	n.a.	n.a.
<b>Negative Predictive Value</b>	n.a.	n.a.
<b>Reference Method</b>	SARS-CoV-2 E-gene assay Thermofisher Scientific	


Victor M. Corman VCH, Tobias Bleicker, Marie Luisa Schmidt, Barbara Mühlemann, Marta Zuchowski, Wendy Karen Jó Lei, Patricia Tscheak, Elisabeth Möncke-Buchner, Marcel A. Müller, Andi Krumbholz, Jan Felix Drexler, Christian Drosten. Comparison of seven commercial SARS-CoV-2 rapid Point-of-Care Antigen tests. medRxiv 2020; medRxiv preprint doi: <https://doi.org/10.1101/2020.11.12.20230292>;

# Corman et al., Germany: Result Details



- a) Distribution of SARS-CoV-2 viral RNA concentrations across clinical samples used for AgPOCT testing.
- b) Overview of tested samples and corresponding outcomes in the seven AgPOCT (per column). Blue fields correspond to a positive AgPOCT result, red fields to a negative result. Empty fields represent samples that were not tested in the corresponding test.

- I: Abbott Panbio™ COVID-19 Ag Rapid Test
- II: RapiGEN BIOCREDIT COVID-19 Ag
- III: Healgen® Coronavirus Ag Rapid Test Cassette (Swab)
- IV Coris Bioconcept Covid.19 Ag Respi-Strip;
- V: Biopharm RIDA®QUICK SARS-CoV-2 Antigen;
- VI NAL von minden; NADAL COVID19-Ag Test;
- VII: Roche/SD Biosensor SARS-CoV Rapid Antigen Test



Victor M. Corman VCH, Tobias Bleicker, Marie Luisa Schmidt, Barbara Mühlemann, Marta Zuchowski, Wendy Karen Jó Lei, Patricia Tscheak, Elisabeth Möncke-Buchner, Marcel A. Müller, Andi Krumbholz, Jan Felix Drexler, Christian Drosten. Comparison of seven commercial SARS-CoV-2 rapid Point-of-Care Antigen tests. medRxiv 2020. medRxiv preprint doi: <https://doi.org/10.1101/2020.11.12.20230292>

# Corman et al., Germany: Summary



## **Aim:**

To provide a reflection of test performance on analytical properties of 7 newly marketed rapid antigen tests during a low SARS-CoV-2 incidence in summer 2020 in the Northern hemisphere

## **Sensitivity:**

Detection range corresponds to ca. 10 million copies per swab and thus corresponds to a concentration that predicts a virus isolation success of ca. 20% in cell culture\*.

## **Hypothesis:**

Taken other data into consideration<sup>1,2,3,4</sup> positive Ag rapid test results indicate large amounts of virus shedding and may thus indicate the time of infectiousness.

\*the numbers are back calculated and inferred from other studies

<sup>1</sup>Wolfel, R et al. Virological assessment of hospitalized patients with COVID-2019. Nature.2020, 581(7809):465-9; <sup>2</sup>van Kampen et al, Shedding of infectious virus in hospitalized patients with coronavirus disease-2019 (COVID-19): duration and key determinants. medRxiv. 2020:2020.06.08.20125310; <sup>3</sup>Perera et al. SARS-CoV-2 Virus Culture and Subgenomic RNA for Respiratory Specimens from patients with mild Coronavirus Disease. Emerg Infect. Dis. 2020;26(11):2701-4. <sup>4</sup>He X et al: Temporal dynamics in viral shedding and transmissibility of COVID-19. Nat. Med. 2020;26(5):672-5

# Mak et al., Hong Kong: Study Summary

## Purpose of the study

- To compare analytical sensitivity and clinical sensitivity for the three commercially available RAD kits.
- Analytical sensitivity for the detection of SARS-CoV-2 virus was determined by limit of detection (LOD) using RT-PCR as a reference method using respiratory specimens from confirmed COVID-19 patients

## Main results

- The LOD of Standard Q was  $10^{-5}$ . The corresponding Ct value for LOD at  $10^{-5}$  was 28.67.
- In the cross-reactivity test using virus isolates, all were tested negative by the RAD kits. Review of the Ct values showed that specimens missed by the RAD kits had relatively high Ct values.

## Specifics

- To determine LOD between different kits, a respiratory specimen was serially diluted and virus concentrations in each dilution were estimated from Ct value
- Specimen: throat saliva, nasopharyngeal swab and throat swab, nasopharyngeal aspirate and different combinations
- Small number of specimen in the subgroups

## Main Conclusions

Although viral culture was not performed in the present study, the Standard Q was 102 fold less sensitive than RT-PCR, it corresponded to the LOD of viral culture based on our results reported previously.

The authors recommended specimens obtained  $\leq 7$  days after symptom onset for use with the Standard Q. Then, the RAD kit can serve as a COVID-19 filter (filtered out of the infected persons and prevent spread to the others).

# Mak et al., Hong Kong: Study Details



## Standard Ag Test

<b>N, PCR + (%)</b>	280 archive specimens (100%)		
<b>Investigated cohort</b>	respiratory specimens from COVID-19 patients collected by the Public Health Laboratory Services Branch (PHLSB) in Hong Kong were retrieved for this evaluation. All of the specimens were confirmed with SARS-CoV-2 infection by RT-PCR as described		
<b>Samples</b>	mainly nasopharyngeal and throat swabs; Samples were mixed in 2 mL of viral transport media (VTM)		
<b>Symptoms</b>	All of the specimens were confirmed with SARS-CoV-2 infection by RT-PCR		
	<b>NP swab &amp; throat swab</b>	<b>NP swab</b>	<b>Throat saliva</b>
Sensitivity overall	71.4 %	65.7%	71,4%
<b>Ct 12.9-18.4</b>	<b>(13-18) 100%</b>	<b>15-18) 100%</b>	<b>(12-18) 100%</b>
<b>Ct 19.8-28.6</b>	<b>(20- 29) 93.8 %</b>	<b>(19-28) 81.3%</b>	<b>(19-29) 88.2%</b>
Ct 29.0-34.2	(29-34) 10%	(29-35) 10%	(29-33) 11.1
<b>Specificity</b>	<b>n.a.</b>		
<b>PPV / NPV</b>	<b>n.a.</b>		
<b>Reference Method</b>	PCR method not clear, most probably in house method, see <a href="https://doi.org/10.1016/j.jcv.2020.104500">https://doi.org/10.1016/j.jcv.2020.104500</a>		

Calculated sensitivity for Ct <29 is 96%

Mak GCK, Lau SSY, Wong KKY, et al. Analytical sensitivity and clinical sensitivity of the three rapid antigen detection kits for detection of SARS-CoV-2 virus. *J Clin Virol.* 2020;133:104684. doi:10.1016/j.jcv.2020.104684

# Mak et al., Hong Kong: Nasopharyngeal Swab



	Standard Q SD Biosensor	Covid-19 Respi Strip Coris	Nadal Covid-19
<b>Sensitivity Ct (mean)</b>	<b>(16.38) 100%</b>	(16.38) 100%	(16.50) 100%
<b>Sensitivity Ct (mean)</b>	<b>(23.44) 81.3 %</b>	(23.44) 31.3 %	(23.31) 56.3 %
<b>Sensitivity CT (mean)</b>	<b>(31.73) 10%</b>	(31.73) 0 %	31.56 0 %
<b>Sensitivity (overall)</b>	<b>65.7%</b>	40 %	51.4 %
<b>Specificity</b>	<b>100%</b>	100%	100%

Mak GCK, Lau SSY, Wong KKY, et al. Analytical sensitivity and clinical sensitivity of the three rapid antigen detection kits for detection of SARS-CoV-2 virus. J Clin Virol. 2020;133:104684. doi:10.1016/j.jcv.2020.104684

# Mak et al., Hong Kong: Nasopharyngeal and Throat Swab



	Standard Q SD Biosensor	Covid-19 Respi Strip Coris	Nadal Covid-19
<b>Sensitivity Ct (mean)</b>	<b>(15.96) 100%</b>	(15.96) 100%	(15.81) 100%
<b>Sensitivity Ct (mean)</b>	<b>(23.72) 93.8%</b>	(23.72) 31.3%	(23.60) 18.8%
<b>Sensitivity CT (mean)</b>	<b>(32.04) 10%</b>	(32.04) 0%	(31.56) 0%
<b>Sensitivity (overall)</b>	<b>71.4%</b>	40%	51.4 %
<b>Specificity</b>	<b>100%</b>	100%	100%

Mak GCK, Lau SSY, Wong KKY, et al. Analytical sensitivity and clinical sensitivity of the three rapid antigen detection kits for detection of SARS-CoV-2 virus. J Clin Virol. 2020;133:104684. doi:10.1016/j.jcv.2020.104684

# Mak et al., Hong Kong: Details



**Table 2**  
Performance characteristics of the rapid antigen detection kits for the presence of SARS-CoV-2 virus in respiratory specimens.

Specimen type <sup>a</sup>	Coris					Nadal					Standard Q				
	Ct value		No. of specimens		sensitivity	Ct value		No. of specimens		sensitivity	Ct value		No. of specimens		sensitivity
	mean	range	tested	positive		mean	range	tested	positive		mean	range	tested	positive	
<b>NPA &amp; TS</b>															
High	16.70	13.90–18.24	9	7	77.8 %	16.27	11.17–18.20	9	9	100 %	16.70	13.90–18.24	9	9	100 %
Normal	24.07	19.90–28.65	16	1	6.3 %	23.77	19.19–28.67	18	3	16.7 %	24.07	19.90–28.65	16	12	75.0 %
Low	31.83	28.69–33.81	10	0	0 %	31.50	29.17–33.69	8	0	0 %	31.83	28.69–33.81	10	0	0 %
All	24.39	13.90–33.81	35	8	22.9 %	23.61	11.17–18.20	35	12	34.3 %	24.39	13.90–33.81	35	21	60.0 %
<b>NPS &amp; TS</b>															
High	15.96	12.94–18.43	9	9	100 %	15.81	11.82–18.49	9	9	100 %	15.96	12.94–18.43	9	9	100 %
Normal	23.72	19.75–28.59	16	5	31.3 %	23.60	19.72–28.29	16 <sup>b</sup>	3	18.8 %	23.72	19.75–28.59	16	15	93.8 %
Low	32.04	29.00–34.24	10	0	0 %	31.96	28.88–33.99	10	0	0 %	32.04	29.00–34.24	10	1	10 %
All	24.10	12.94–34.24	35	14	40 %	23.98	11.82–33.99	35	12	34.3 %	24.10	12.94–34.24	35	25	71.4 %
<b>NPS</b>															
High	16.38	14.56–18.07	9	9	100 %	16.05	13.30–17.84	9	9	100 %	16.38	14.56–18.07	9	9	100 %
Normal	23.44	18.98–27.51	16	5	31.3 %	23.31	18.89–28.38	16	9	56.3 %	23.44	18.98–27.51	16	13	81.3 %
Low	31.73	28.91–34.47	10	0	0 %	31.56	28.81–33.80	10	0	0 %	31.73	28.91–34.47	10	1	10 %
All	23.99	14.56–34.47	35	14	40 %	23.80	13.30–33.80	35	18	51.4 %	23.99	14.56–34.47	35	23	65.7 %
<b>throat saliva</b>															
High	15.30	11.46–18.37	9 <sup>b</sup>	8	88.9 %	15.27	11.32–18.36	9 <sup>c</sup>	7	77.8 %	15.30	11.46–18.37	9	9	100 %
Normal	23.21	18.68–28.62	17 <sup>b</sup>	3	17.6 %	23.20	18.66–28.62	17	4	23.5 %	23.21	18.68–28.62	17 <sup>d</sup>	15	88.2 %
Low	31.06	29.18–33.31	9	0	0 %	31.04	29.16–33.30	9 <sup>b</sup>	0	0 %	31.06	29.18–33.31	9	1	11.1 %
All	23.20	11.46–33.31	35	11	31.4 %	23.18	11.32–33.30	35	11	31.4 %	23.20	11.46–33.31	35	25	71.4 %

<sup>a</sup> ‘High’, means specimens with Ct values <18.57 of SARS-CoV-2 virus RT-PCR; ‘Normal’, Ct values between 18.57 and 28.67; ‘Low’, Ct values >28.67.

<sup>b</sup> The control band was failed to appear for one specimen.

<sup>c</sup> The control band was failed to appear for two specimens.

<sup>d</sup> The nozzle cap was blocked for one specimen, 100 µL of suspension was added directly to the device. The final test result was positive.

# Chaimayo et al., Thailand: Study Summary



## Purpose of the study

Performance characteristics of the rapid SARS-CoV-2 antigen test were evaluated and compared with the gold standard RT-PCR for diagnosis of COVID-19 cases.

## Main results

The rapid assay for SARS-CoV-2 antigen detection showed comparable sensitivity and specificity with the RT-PCR assay.

- Sensitivity 98.33%
- Specificity 98.73%

## Specifics

Cohort: suspected COVID-19 cases, including pre-operative patients. Mainly combined nasopharyngeal and throat swabs were used.

## Main Conclusions

The rapid SARS-CoV-2 antigen test can benefit all healthcare workers in managing infected individuals in time effectively, in high prevalence areas and especially in rural and outbreak areas. The advantage of the Standard Q COVID-19 Ag test as a screening for COVID-19 is its simple procedure and quick results with high NPV, but its disadvantage is low PPV in a low prevalence area.

# Chaimayo et al, Thailand: Study Details



## Standard Ag Test

<b>N, PCR + (%)</b>	454 <b>(13.2%)</b>
<b>Investigated cohort</b>	suspected COVID-19 cases, including pre-operative patients
<b>Samples</b>	mainly nasopharyngeal and throat swabs; Samples were mixed in 2 mL of viral transport media (VTM)
<b>Symptoms</b>	three days (range 0–14),
<b>Sensitivity overall</b>	<b>98.33%</b> (95% CI, 91.06–99.96%) One negative sample had Ct values of E, RdRp, and N with 31.08 / 39.2 / 35.54 (negative RT-PCR is defined as having Ct-values larger than 40)
<b>Specificity</b>	<b>98.73%</b> (95% CI, 97.06–99.59%)
<b>PPV / NPV</b>	PPV and NPV of the assay could not be accurately calculated without the present population prevalence of COVID-19.
<b>Reference Method</b>	Allplex™ 2019-nCoV Assay (Seegene®, Korea)

## Purpose of the study

A manufacturer-independent, prospective diagnostic accuracy study with comparison of a supervised, self-collected anterior nose (AN) swab sample with a professional collected nasopharyngeal swab (NP) sample, using STANDARD Q COVID-19 Ag Test (SD Biosensor)

## Main results

The Ag-RDT with AN sampling showed a sensitivity of 74.4% and specificity of 99.2% compared to RT-PCR. The sensitivity with NP sampling was 79.5% and specificity was 99.6%. In patients with high viral load ( $>7.0$  log<sub>10</sub> RNA SARSCoV2/swab), the sensitivity of the Ag-RDT with AN sampling was 96% and 100% with NP sampling.

## Specifics

A supervised self-collected nasal sample (both nostrils) were taken first, then the combined NP/OP (1 nostril) for PCR, lastly the NP (the other nostril) for the Ag test was taken. Sequence might lead to different viral loads. NP swab was usually rotated against the nasopharyngeal wall for **less** time than recommended by the manufacturer

## Main Conclusions

- Supervised self-sampling from the anterior nose is a reliable alternative to professional nasopharyngeal sampling using a WHO-listed SARS-CoV-2 Ag-RDT
- The Ag-RDT frequently did not detect patients with lower viral load or with symptoms  $>7$  days

# Lindner et al., Germany: Study Details



## Roche Rapid Ag Test

**N, PCR + (%)**

289 (13.5%)

**Investigated cohort**

Adults at high risk according to clinical suspicion  
On the day of testing, 97.6% of participants had one or more symptoms consistent with COVID-19.

**Samples**

Supervised anterior nose swab (AN) -- > off-label

**Professional NP swab**

**Symptoms**

Average 4.4 days (SD 2.7)

**Sensitivity overall**

**74.4% (CI 58.9-85.4)**

**79.5 (CI 64.5-89.2)**

**Sensitivity high viral load  
(>7.0 log<sub>10</sub> RNA SARS-CoV2/swab)**

96% (CI 80.5-99.3)

**100% (CI 86.7-100)**

**Ct 17.3-23.7**

95.7%

**100%**

**Ct 17.3-25.3**

92.3 %

**96.2%**

**Ct 17.3-29.6**

87.1%

**90.3 %**

**Ct 17.3-30.0**

84.4%

**87.5%**

**Ct 24.2-35.5**

43.8%

**50.0%**

**Ct 25.3- 35.5**

38.5%

**46.2%**

Sensitivities calculated based on Table in the publication

**Specificity**

99.2% (CI 97.1-99.8)

**99.6 (CI 97.8-100)**

**Pos % agreement AN / NP**

**90.6% (Ci 75.8-96.8)**

**Reference Method**

The Roche cobas SARS-CoV-2 assay or the SARS-CoV-2 E-gene assay from TibMolbiol (Berlin, Germany)

Lindner et al 2020 doi: <https://doi.org/10.1101/2020.10.26.20219600>

# Igloi et al., The Netherlands: Study Summary



## Purpose of the study

The Roche/SD Biosensor lateral flow antigen rapid test was evaluated in a mild symptomatic population at a large drive through testing site.

## Main results

Overall sensitivity and specificity were 84.9% and 99.5%  
Sensitivity for samples with high loads of viral RNA (ct <30, 2.17E+05 E gene copy/ml) and who presented within 7 days since symptom onset increased to 95.8% .

## Specifics

All Ag Rapid Antigen Tests and PCR positive samples were cultured to correlate results with infectivity. Eligibility for a free of charge test includes either symptoms or close contact with a confirmed SARS-CoV-2 infected person, therefore symptoms may be over-reported.

## Main Conclusions

- People with early onset and high viral load were detected with 98.2% sensitivity, 97% of individuals in which virus could be cultured were detected by the rapid test.
- This test is suitable to detect mild symptomatic cases, suggesting screening based on Ag RDT alone in this population would have a high sensitivity for ruling out infectious individuals .

# Igloi et al., The Netherlands: Study Details



## Roche Rapid Ag Test

<b>N, PCR + (%)</b>	970 (19.2%)		
<b>Investigated cohort</b>	Mild symptomatic population, eligibility for a free of charge test includes either symptoms or close contact with a confirmed SARS-CoV-2 infected person		
<b>Samples</b>	First swab: combined NP + OP for PCR and viral cell culture; in UTM (HiViral™) Nasopharyngeal swabs for Rapid Ag Test as a second swab from the same nostril		
<b>Symptomatics, n (%)</b>	(xx%)		
DPSO (median)	4		
Days < 0-3)	44.0%		
Days 4-7	45.7%		
Days 8+	10.3%		
<b>PCR Ct (median; CI)</b>	23.6 (15.6-37.4)		
	<b>0-3 days post onset</b>	<b>0-7 days post onset</b>	<b>All</b>
Clinical Sensitivity	94.9 (86.1-98.3), 319	90.6 (84.3-94.6), 650	84.9 (79.1-89.4), 970
Sensitivity CT < 30 (95% CI), N	98.2 (90.6-99.9), 316	95.8 (90.5-98.2), 640	94.3 (89.6-99.9), 943
Sensitivity CT < 25 (95% CI)	100 (92.1-100), 305	98.8 (93.7-99.9), 608	99.1 (95.2-100), 897
<b>PPV</b>	98.2 (90.7-99.9)	98.3 (94.0-99.5)	97.5 (93.8-99.0)
<b>Clinical specificity (95% CI), N</b>	99.6 (97.9-100), 319	99.6 (98.6-99.9), 650	99.5 (98.7-99.8), 970
<b>Reference Method</b>	cobas 6800 and Vero cell clone 118; sample material: combined NP + OP swabs		

# Krüttgen et al., Germany: Study Summary



## Purpose of the study

The sensitivity and specificity of the new Roche SARS-CoV-2 Rapid Antigen Test was evaluated

## Main results

- The assay's sensitivity with samples with a cycle threshold of  $< 25$  was 100% and gradually decreases to 22,2% with cycle thresholds  $\geq 35$ .
- They found a specificity of 96%.
- Samples with Ct-values  $>30$  usually do not allow culturing of the virus indicating low infectivity.

## Specifics

Using 75 swabs from patients previously tested positive by SARS-CoV-2 PCR and 75 swabs from patients previously tested negative by SARS-CoV-2 PCR,

## Main Conclusions

Sensitivity and specificity of the antigen assay is inferior to the PCR assay, but the overall sensitivity is strictly dependent on the distribution of cycle thresholds (Ct) within the population of specimens and does not allow a realistic evaluation of the assay. The new test might be useful to rapidly identify contagious individuals as they state that samples with Ct-values  $>30$  usually do not allow culturing of the virus indicating low infectivity.

Krüttgen A, Cornelissen CG, Dreher M, Hornef MW, Im"ohl M, Kleines M, Comparison of the SARS-CoV-2 Rapid Antigen Test to the Real Star Sars-CoV-2 RT PCR Kit, *Journal of Virological Methods* (2020), doi: <https://doi.org/10.1016/j.jviromet.2020.114024>

# Krüttgen et al., Germany: Study Details



## Roche Ag Test

<b>N, PCR + (%)</b>	<b>150 (50%)</b> (selected samples)
<b>Investigated cohort</b>	Using 75 swabs from patients previously tested positive by SARS-CoV-2 PCR and 75 swabs from patients previously tested negative by SARS-CoV-2 PCR
<b>Samples</b>	350 µl of swab transport medium were mixed with extraction buffer provided by the manufacturer
<b>Symptoms</b>	n.a.; sample collection contained clinical specimens only and the SARS-CoV-2 RNA positive subpopulation was characterized by a wide range of Ct-values with medium and low Ct-values dominating.
<b>Sensitivity overall</b>	70,7%
<b>Sensitivity Ct &lt; 20</b>	100%
<b>Sensitivity Ct 25-30</b>	95%
<b>Sensitivity Ct 30-35</b>	44.8%
<b>Sensitivity Ct &gt;35</b>	22.2%,
<b>Specificity</b>	<b>96%</b>
<b>Reference Method</b>	Real Star SARS-CoV-2 RT PCR Kit (Altona, Germany)

Krüttgen A, Cornelissen CG, Dreher M, Hornef MW, Imohl M, Kleines M, Comparison of the SARS-CoV-2 Rapid Antigen Test to the Real Star Sars-CoV-2 RT PCR Kit, *Journal of Virological Methods* (2020), doi: <https://doi.org/10.1016/j.jviromet.2020.114024>

# Nalumansi et al., Uganda: Study Summary



## Purpose of the study

- The aim of this study was to evaluate a low cost, easy-to-use rapid antigen test for diagnosing COVID-19 at the point-of-care.
- Ag Test and results compared with the qRT-PCR results

## Main results

- Sensitivity and specificity of the antigen test were 70.0% (95% CI: 60 - 79) and 92% (95% CI: 87 - 96) respectively; diagnostic accuracy was 84% (95% CI: 79 - 88).
- The antigen test was more likely to be positive in samples with qRT-PCR Ct values  $\leq 29$  reaching a sensitivity of 92%.

## Specifics

- Nasopharyngeal swabs from suspect COVID-19 cases and from low-risk volunteers were tested on the STANDARD Q COVID-19
- 262 samples incl 90 RT-PCR positives
- The sequence of sampling is not clear

## Main Conclusions

- They conclude that the STANDARD Q COVID-19 Ag Test performed less than optimally in this evaluation but that it may still have an important role to play early in infection when timely access to molecular testing is not available but results should be confirmed by qRT-PCR.
- “Unusual” categorization of the Ct values: they were categorized as strongly positive (Ct  $\leq 29$ ) (indicative of abundant target nucleic acid in the sample), moderately positive (Ct 30-37) and weakly positive (Ct 38-39)

<https://doi.org/10.1016/j.ijid.2020.10.073> IJID 4794

# Nalumansi et al., Uganda: Study Details



## Roche Ag Test

<b>N, PCR + (%)</b>	<b>262 (34.4%)</b>
<b>Investigated cohort</b>	suspect COVID-19 cases and from low-risk volunteers were tested on the STANDARD Q COVID-19, 262 samples incl. 90 RT-PCR positives
<b>Samples</b>	Nasopharyngeal swabs
<b>Symptoms</b>	n.a., 14% of the positives were mildly symptomatic – no data on symptom onset
<b>Sensitivity overall Ct ≤29-39</b>	70% (95% CI: 60 - 79)
<b>Sensitivity Ct ≤29</b>	92% (95% CI: 87- 96)
<b>Sensitivity Ct 30-37</b>	55%
<b>Sensitivity Ct 38-39</b>	56%
<b>Specificity</b>	<b>92%</b> (95%CI 87-96)
<b>Reference Method</b>	Berlin protocol for RT-PCR

Ct values not well comparable with other studies

# Schwob et al., Switzerland: Study Summary



## Purpose of the study

A prospective clinical trial in symptomatic patients to investigate analytical (PCR and RDTs) and sampling procedures (saliva and NP swab) and in order to compare the detection rate of SARS-CoV-2 and sensitivities of i) RDT on NP swab, ii) PCR on NP swab and iii) PCR on saliva.

Secondary objectives were to compare detection rates and sensitivities stratified by Viral Load (VL) categories.

## Main results

The results of the present study show that the detection rate of positive COVID-19 cases by RDT was high, especially for those with a VL of  $\geq 10^6$  copies/ml.

There was a slight variability in performance between the three different RDTs with STANDARD Q® having a higher sensitivity (93%) than those of Panbio™ (86%) and COVID-VIRO® (84%).

## Specifics

Very low inter-observer variation in test line reading which confirms user-friendliness.

Well defined population presenting within 7 days after symptom onset.

## Main Conclusions

The high performance of RDTs allows rapid identification of COVID cases with immediate isolation of the vast majority of contagious individuals. Based on the 100% specificity of high quality RDT there is no need to confirm a positive RDT test result by an additional PCR test.

# Schwob et al., Switzerland: Study Details



## Roche SARS-CoV-2 Rapid Ag

## Panbio Abbott

## Coivid-Viro Ag tests

**N, PCR positive (%)**

928 (40.1% (36.9–43.3%) by NP PCR)

**Investigated cohorts**

96% of participants had at least one major symptom and 4% at least one minor and a close contact with a documented COVID-19 case. Mean duration of symptoms at the time of swab collection/testing was 2.6 days (SD 2.3, range 0-30).

**Samples**

two nasopharyngeal swabs, one for PCR and one for RDT analyses (sequence not described)

**Sensitivity (95% CI)**

**92.9%** (86.4–96.9)

86.1% (78.6–91.7%)

84.1% (76.9–89.7%)

Ct ≤26 or VL\* ≥ 10<sup>6</sup> (Ct<sup>26</sup>), (95%CI)

96.6% (90.5–99.3)

97.8% (92.1–99.7%)

95.3% (89.4–98.5%)

**Specificity (95%CI)**

100% (99.3–100)

**Reference Method**

in-house RT-PCR on the automated molecular diagnostic platform targeting the E gene,13–15 or using the SARS-CoV-2 test of the Cobas 6800 instrument (Roche, Basel, Switzerland).

\*The thresholds chosen for analyses by VL were 10<sup>5</sup> copies/ml (Ct=30) and 10<sup>6</sup> copies/ml (Ct=26), based on recent and older data investigating the link between viral loads and the presence of culture-competent virus<sup>1-5</sup>

# Salvagno et al., Italy: Study Summary



## Purpose of the study

The purpose of this study was the clinical assessment of the new Roche SARS-CoV-2 Rapid Antigen Test versus a PCR assay in nasopharyngeal swabs.

## Main results

The sensitivity was found to range between 97-100% in clinical samples with Ct values <25, between 50-81% in those with Ct values between 25-<30, but low as 12-18% in samples with Ct values between 30-<37.

## Specifics

The study population consisted of all consecutive patients referred for SARS CoV- 2 diagnostic testing to the Hospital.

## Main Conclusions

The clinical performance of Roche SARS-CoV-2 Rapid Antigen Test is excellent in nasopharyngeal swabs with Ct values <25, which makes it a reliable screening test in patients with high viral load.

# Salvagno et al., Italy: Study Details



## Roche Ag Test

<b>N, PCR + (%)</b>	<b>321 (46.4%)</b>
<b>Investigated cohort</b>	The study population consisted of all consecutive patients referred for SARS CoV-2 diagnostic testing to the Pederzoli Hospital;
<b>Samples</b>	A single swab (Virus swab UTM™, Copan, Brescia, Italy) was collected from each patient and concomitantly used for both Roche SARS-CoV-2 Rapid Antigen testing and molecular testing in 350 µl volume.
<b>Symptoms</b>	n.a.
<b>Sensitivity overall</b>	72.5%
<b>Sensitivity Ct &lt; 25</b>	97-100%
<b>Sensitivity Ct 25-&lt;30</b>	50-81%
<b>Sensitivity Ct 30-37</b>	12-18%
<b>Specificity</b>	<b>99.4%</b>
<b>Reference Method</b>	Seegene Allplex™2019-nCoV Assay, (Seegene, South Korea), targeting three viral genes (N, E and RdRP),

Salvagno GL, Gianfilippi G, Bragantini D, Henry BM, Lippi G. Clinical assessment of the Roche SARS-CoV-2 Rapid Antigen Test. *Diagnosis (Berl)*. 2020. doi: 10.1515/dx-2020-0154

# External Clinical Performance Study Results Overview

## *Roche SARS-CoV-2 Rapid Antigen Test*

Study	#Sample	# PCR+ (%)	Sensitivity (95% CI) Ct ≤ x	Sensitivity (CI)	Specificity (CI)
<b>FIND, BRA &amp; D</b>	1659	9.2%	<b>97.14%</b> (90.1-99.65) <b>Ct≤25</b>	84.97% (78.3-90.23)	<b>98.94%</b> (98.23-99.39)
<b>HUG, CH</b>	529	36%	<b>98%</b> (n.a.) <b>Ct≤22</b>	89.0% (83.69-93.06)	<b>99.70%</b> (98.36-99.99)
<b>Cerutti, I</b>	330	33%	<b>100%</b> (n.a.) <b>Ct≤28</b>	72.1% (83.69-93.06)	<b>100%</b> (98.36-100)
<b>Krueger, D &amp; UK</b>	1263	3%	<b>100%</b> (82.4-100) <b>Ct≤25</b>	76.6% (62.8-86.4)	<b>99.3%</b> (98.6-99.6)
<b>Van Beek, NL</b>	1754	100%	Detection of culture positive and RT-PCR-confirmed: <b>94.3-99.8%</b>		
<b>Corman, D</b>	115	n.a.	6.78 copies/swab LoD, 95% mean hit rate detected as little as 4.4 PFU (plaque forming units) of virus per test.		<b>97.12%</b> n= 35 Cumulative Spec. 98.53%
<b>Mak, HK</b>	280	100%	<b>96%</b> <b>Ct&lt;29</b>	<b>71.4%</b>	<b>n.a.</b>

**CT-values cannot be compared 1:1 as RT-PCR methods vary across sites with different genome targets, PCR instruments and reagents**

# External Clinical Performance Study Results Overview

## *Roche SARS-CoV-2 Rapid Antigen Test*

Study	#Sample	# PCR+ (%)	Sensitivity (95% CI) Ct ≤ x	Sensitivity (95% CI)	Specificity (95% CI)
<b>Chaimayo, TAI</b>	454	13.2%	<b>98.3%</b> (91.06–99.96%) Ct n.a.	98.3% (95% CI, 91.06–99.96%)	98.7% (97.06–99.59%)
<b>Lindner, D</b>	289	13.5%	<b>96.2%</b> Ct <b>17.3-25.3</b>	74.4% (CI 58.9-85.4)	99.6 (CI 97.8-100)
<b>Igloi, NL</b>	970	19.2%	<b>99.1%</b> (95.2-100) Ct < <b>25</b>	84.9 (79.1-89.4)	99.5 (98.7-99.8)
<b>Krüttgen, D</b>	150	50%	<b>100%</b> Ct < <b>25</b>	70.7%	96%
<b>Nalumansi; UG</b>	262	34.4%	<b>92%</b> Ct ≤ <b>29</b>	70%	92% (95%CI 87-96)
<b>Schwob, CH</b>	928	40.1%	<b>96.6%</b> (90.5-99.3) Ct ≤ <b>26</b>	92.9% (86.4-96.9)	100%
<b>Salvagno, I</b>	321	46.4%	<b>97-100%</b> Ct < <b>25</b>	72.5%	99.4%

**CT-values cannot be compared 1:1 as RT-PCR methods vary across sites with different genome targets, PCR instruments and reagents**

# Conclusions



- 14 studies presented with over 9'300 patient samples
- The sensitivity of the Roche / SD Biosensor POC Antigen assay was between 96.2 to 100% with a CT that is considered to be associated with culture positive results. \*
- If the specimens are obtained  $\leq 7$  days after symptom onset for use with the Rapid Antigen test, it can help to filter out the infected persons and prevent spread to the others.
- First real world performance data confirms the primary use case for POC assay, however, more and larger studies are needed.

\*The data from Uganda are not considered due to great discrepancy of the Ct values and categorization compared to all other republications.

# Thank you

*for your attention*

---

COBAS is a trademark of Roche.

© 2020 Roche

**Published by:**

Roche Diagnostics International Ltd

CH-6343 Rotkreuz

Switzerland

**[diagnostics.roche.com](https://diagnostics.roche.com)**

*This presentation is our intellectual property. Without our written consent, it shall neither be copied in any manner, nor used for manufacturing, nor communicated to third parties.*

# *Doing now what patients need next*